



Glass Availability Guide

Single Vision		
Single Vision	Color	Base Curve
66/65	White	Plano, 2.25, 4.25, 6.25, 8.25, 10.25
	PGX	
	PBX	4.25, 6.25, 8.25
	White 1.60 Index (70mm)	Plano, 2.25, 4.25, 6.25
	White 1.70 Index	Plano, 2.25, 4.25
70/71mm	White	Plano, 2.25, 4.25, 6.25, 8.25
	PGX	
	Thin & Dark	
70mm	Gray 3 Polar	4.00, 5.00, 6.00, 8.00
	Brown 3 Polar	
71mm	White	Plano, 2.25, 4.25, 6.25, 8.25
	PGX	
	PBX	
	Thin & Dark	
	Grey C Polar	3.25, 4.25, 6.25, 8.25
	Brown C Polar	
	G-15	4.25, 6.25, 8.25
	1.8 White	Plano, 2.25, 4.25
	X-Ray 1.80 Filter .75pB	4.00, 6.00, 8.00
	76mm	White
PGX		

Bifocals

Bifocals				
D28	Color	Base Curve	D28	
65/66mm	White	2.25, 4.25, 6.25, 8.25	Base	ADD range
	PGX		plano	.75 - 3.75
	PBX			
70/71mm	White	Plano, 2.25, 4.25, 6.25 8.25, 10.25	2.25	.75 - 4.50
	PGX		4.25	.50 - 5.00
	PBX		6.25	.50 - 5.75
	Thin & Dark			
	X-Ray 1.80 Filter .75pB		4.00, 6.00, 8.00	
75mm	White	4.25, 6.25	8.25	.50 - 6.50
	PGX			
76mm	White	Plano, 2.25	10.25	.50 - 7.00
	White	4.25, 6.25, 8.25		
	PGX			

D35	Color	Base Curve	D35	
66mm	White	4.25, 6.25, 8.25	Base	ADD range
	PGX		plano	.75 - 2.75
70/71mm	White	4.25, 6.25, 8.25	2.25	.75 - 3.50
	PGX		4.25	.50 - 4.00
71mm		4.25, 6.25, 8.25	2.25	.50 - 3.50
			6.25	.50 - 4.25
	White			
	PGX			
	Fixed Tint			
76mm	White	4.25, 6.25, 8.25	8.25	.50 - 5.50
	PGX		10.25	.50 - 6.00

R22	Color	Base Curve	R22	
66mm	White	4.25, 6.25, 8.25	Dia	ADD range
	PGX		66	1.00 - 3.00
71mm	White	Plano, 2.25, 4.25, 6.25, 8.25	71	.50 - 3.00
	PGX			

Trifocals and Occupational

Trifocals and Occupational				
7x28	Color	Base Curve	7x28	
65/66mm	White	4.25, 6.25, 8.25	<i>Base</i>	<i>ADD range</i>
	PGX		plano	1.50 - 3.75
70/71mm	White	4.25, 6.25, 8.25 / Plano, 2.25, 10.25	2.25	1.50 - 4.50
	PGX		4.25	1.25 - 5.00
71mm	White	Plano, 2.25, 4.25, 6.25, 8.25, 10.25	6.25	1.00 - 5.00
	PGX		8.25	1.25 - 6.50
	PBX		10.25	1.50 - 5.00
	Thin & Dark		Plano, 2.25, 4.25, 6.25, 8.25	
76mm	White	4.25, 6.25, 8.25		
	PGX			

7x35	Color	Base Curve	7x25 / 7x28	
71mm	White	4.25, 6.25, 8.25	<i>Base</i>	<i>ADD range</i>
	PGX		ALL	1.00 - 3.00

DD28	Color	Base Curve	DD28	
71mm	White	4.25, 6.25, 8.25	<i>Base</i>	<i>ADD range</i>
	PGX		4.25	1.00 - 3.00
			6.25	1.00 - 3.00
76mm	White	4.25, 6.25, 8.25	8.25	1.00 - 3.00
	PGX			

DD35	Color	Base Curve	DD35	
			<i>Base</i>	<i>ADD range</i>
76mm	White	4.25, 6.25, 8.25	All	1.00 - 3.00
	PGX			

D & 7x28	Color	Base Curve	Q28	
71mm	White	4.25, 6.25, 8.25 (limited availability)	<i>Base</i>	<i>ADD range</i>
	PGX		All	1.00 - 3.00

7x35	Color	Base Curve		
76mm	PGX	4.25, 6.25, 8.25	<i>Base</i>	<i>ADD range</i>
			All	1.50 - 3.00

Progressives				
Freedom ID Progressive	Color	Base Curve	Freedom ID	
71mm	White	3.25, 5.25, 7.25	<i>Base</i>	<i>ADD range</i>
	PGX		All	1.00 - 3.50
	PBX			
	Thin & Dark			
	Gray 3 Polar			
	Brown 3 Polar			
Essilor Adaptar	Color	Base Curve	Adaptar	
70mm	White	2.25, 4.25, 6.25, 8.25	<i>Base</i>	<i>ADD range</i>
	PGX		All	1.00 - 3.50
Essilor Natural 1.60	Color	Base Curve	Natural 1.60	
70mm	White	3.25, 5.25, 6.25	<i>Base</i>	<i>ADD range</i>
	PGX		All	1.00 - 3.50
Varilux Comfort 1.60	Color	Base Curve	Varilux Comfort 1.60	
70mm	White	3.25, 4.25, 5.25, 7.25	<i>Base</i>	<i>ADD range</i>
	PGX		All	1.00 - 3.50
Sola VIP	Color	Base Curve	Sola VIP	
70mm	White	2.25, 3.25, 4.25, 5.25 7.25-65mm only	<i>Base</i>	<i>ADD range</i>
	PGX		All	1.00 - 3.50
Zeiss GT2	Color	Base Curve	Zeiss GT2	
71mm	White	3.25, 4.25, 5.25 6.25-65mm Only	<i>Base</i>	<i>ADD range</i>
	PGX		All	1.00 - 3.50
Digital Progressive	Color	Base Curve		
71mm	X-Ray 1.80 Filter .75pB	4.00, 6.00, 8.00		