

Discover the future of glass

crystalCLEAR digital glass from Vision Dynamics Lab

The benefits of digital are undeniable

- Advanced prescription accuracy
- Exceptional optical clarity
- Optimized center thickness
- Suitable for the majority of prescriptions

Combined with the benefits of

- High-performance glass for a thinner, lighter lens
- Superior scratch resistance and chemical durability
- Can be surfaced 25% thinner than traditional glass
- Proprietary molecular structure



crystalCLEAR digital glass

AVAILABILITY	crystalCLEAR Single Vision -7.00 to +5.00 cyl. to -4.00	crystalCLEAR Progressive -7.00 to +5.00 cyl. to -4.00 +0.25 to +4.50 add
1.523 Clear	✓	✓
1.523 Thin & Dark	✓	✓
Edging & Ion-Exchange Processing	Included	Included
ADD-ONS		
Sphere Powers > +/-4.00 (per lens per diopter)		Add +3.25 to +4.00
Cylinder Powers > -3.00 (per lens per diopter)		Add +4.25 to +4.50
Lenticularization		Prism (per lens per diopter)
AR Coating		Handling
		Shipping: Actual or UPS Account #

crystalCLEAR Progressive
General use progressive design with a smooth transition between distance and near vision offering the easiest way to find focus points. The crystalCLEAR progressive has high performance and easy adaptation.

Minimum fitting height 14mm